


MEDICA EASYCAL DRUGS OF ABUSE



REF 14696 1 X 5 ML BOTTLE, CALIBRATOR

ETHYL ALCOHOL (ETOH) CALIBRATOR

The bottle contains 5.0 mL aqueous based ready-to-use liquid spiked with a known concentration of Ethyl Alcohol (Ethanol).

INTENDED USE

The Medica Ethyl Alcohol Calibrator is intended for in-vitro diagnostic use for the calibration of the ethyl alcohol assay, which is used on the EasyRA for the quantitative determination of ethyl alcohol in urine.

SAFETY PRECAUTIONS AND WARNINGS

1. The human based material from which this product has been derived has not been tested for the presence of antibodies to HIV, the hepatitis antigens, or anti-hepatitis antibodies. This material and all patient samples should be handled and disposed of as though they are capable of transmitting infectious diseases.
2. The Ethyl Alcohol Calibrator contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing, always flush with copious volumes of water to prevent azide build-up
3. Do not use the calibrator beyond the expiration date.
4. The bottle containing this calibrator should be tightly capped and refrigerated when not in use to prevent evaporation.

INSTRUCTIONS FOR HANDLING, STORAGE AND STABILITY

The calibrator material is stable until the expiration date on the label if stored at 2° - 8°C. The calibrator should not be frozen or exposed to temperatures above 32°C.

INSTRUCTIONS FOR USE

The calibrator material is ready-to-use. Remove bottle cap. Dispense three drops of material into a sample cup and proceed with instructions in the EasyRA Operator's Manual.

VALUE ASSIGNMENT

The table below contains the assigned values for the Ethyl Alcohol Calibrator:

Lot No. XXXXX Expires: XXXX-XX-XX

Analyte:	Assigned Value	Units
Ethyl Alcohol	100	mg/dL