

MEDICA EASYCAL DRUGS OF ABUSE

REF 14799 1 X 5 ML BOTTLE, CALIBRATOR

UNIVERSAL NEGATIVE CALIBRATOR, 0 ng/mL

The bottle contains human urine-based, drug-free, ready-to-use liquid.

INTENDED USE

The Medica Universal Negative Calibrator, in combination with the Cutoff Calibrators, is intended for in-vitro diagnostic use for the calibration of drugs of abuse assays, which is used on the EasyRA for the qualitative determination of specific drugs of abuse.

SAFETY PRECAUTIONS AND WARNINGS

1. The human based material from which this product has been derived has not been tested for the presence of antibodies to HIV, the hepatitis antigens, or anti-hepatitis antibodies. This material and all patient samples should be handled and disposed of as though they are capable of transmitting infectious diseases.
2. The calibrator contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing, always flush with copious volumes of water to prevent azide build-up.
3. Do not use the calibrator beyond the expiration date.
4. The bottle containing this calibrator should be tightly capped and refrigerated when not in use to prevent evaporation.

INSTRUCTIONS FOR HANDLING, STORAGE AND STABILITY

The calibrator material is stable until the expiration date on the label if stored at 2- 8°C. The calibrator should not be frozen or exposed to temperatures above 32°C.

INSTRUCTIONS FOR USE

The calibrator material is ready-to-use. The Medica EasyCal Universal Negative Calibrator, REF 14799 and a Medica Cutoff Calibrator are required for the calibration of the assay. Refer to assay and cutoff calibrator instructions for use. The calibrators should be placed in 2 different positions on the sample wheel in ascending concentration. Remove bottle cap. Dispense three drops of material into a sample cup and proceed with instructions in the EasyRA Operator's Manual.

VALUE ASSIGNMENT

The table below contains the assigned values for the Universal Negative Calibrator:

Lot No. XXXXX Expires: YYYY-MM-DD

Assigned Value	Units
----------------	-------

0	ng/mL
---	-------

Manufactured for: