

REF 10764

1 x 29 mL/10 mL

Precision Test Dye / Präzisionstest-Farbstoff / Test de precisión de colorante / Colorant pour test de précision / Hassaslık testi boyası / Colorante per test di precisione / Corante para teste de precisão / Краситель для анализа точности / Barwnik do testu precyzji / 精密度试验用染液

Precision Test Dye
Each dual reagent wedge contains two levels of the dye solution.
The large section of the wedge contains 29 mL; the smaller section of the wedge contains 10 mL.

Intended Use
For use in testing the precision of pipetting sample at 2 µL and 16 µL for the Medica EasyRA Chemistry Analyzer and as a diagnostic tool.
The dye in the large section of the wedge is more dilute and is used to test the precision of a 16 µL sample. The dye in the small section of the wedge is used to test the precision of a 2 µL sample.

Precautions
Use gloves when handling and dispose of in a safe manner. Avoid eye and skin contact. Wash hands after using this material.

Reagents
Large wedge section: potassium dichromate at 10 mmol/L in 0.01 N Sulfuric acid.
Small wedge section: potassium dichromate at 100 mmol/L in 0.01 N Sulfuric Acid.

Instructions For Handling, Storage And Stability
The precision test solutions are ready-to-use as supplied.
The precision dye wedge is stable to the expiration date on the label if stored at room temperature (18°–25°C) in the closed wedge.
The solutions should be clear and yellow (large side), orange (small side). Do not use if any turbidity or precipitate appears as this may indicate contamination.

Instructions for Use
It is recommended that precision testing be performed on a weekly basis and can be used as a diagnostic troubleshooting test.
NOTE: Gently invert the wedge three or four times before placing in the reagent tray.
Remove both caps and place the Precision wedge in the EasyRA reagent tray. Verify that there are at least 20 unused cuvettes in the reagent area.
Refer to the EasyRA Operators Manual for complete instructions on performing the Precision Test.

Expected Values
The precision dye solution at 100 mmol/L is expected to have an absorbance in the range of 1.1-1.5 when 2 µL of dye is mixed in a cuvette with 255 µL of diluent and analyzed at 340 nm. The CV specification is < 2%.
The precision dye solution at 10 mmol/L is expected to have an absorbance in the range of 1.1-1.4 when 16 µL of dye is mixed in a cuvette with 180 µL of diluent and analyzed at 340 nm. The CV specification is < 1%.

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Medica Corporation, 5 Oak Park Drive